ELTEN ESD Pro
HAIX® CO SYSTEM
ATLAS KLIMA-KOMFORT®
STEITZ SECURA ERGO-SOFT ESD
FERTIGUNGSANWEISUNG

FACHHANDEL

EN – ELTEN ESD Pro/HAIX® CO SYSTEM/ATLAS KLIMA-KOMFORT®/

Steitz Secura ERGO-SOFT ESD –

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T. +49 (0)7021 738 30-0
info@perpedes.de
www.perpedes.de

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CONTENTS

Introduction	22
Production instruction ELTEN ESD Pro/ Type-approved treatment system and its intended purpose..	23
Indications/contraindications/processing instructions.....	24
Assembly and adaption regulations/ labeling of the custom-made device	25
Production instruction HAIX® CO system/ Type-approved treatment system and its intended purpose..	26
Indications/contraindications/processing instructions.....	27
Assembly and adaption regulations	28
Labeling of the custom-made device	29

For the visually impaired, this document is available in PDF format with this QR code.

Production instruction ATLAS KLIMA-KOMFORT®/ Type-approved treatment system and its intended purpose..	38
Indications/contraindications/processing instructions.....	39
Assembly and adaption regulations	40
Labeling of the custom-made device	41
Production instruction Steitz Secura ERGO-SOFT ESD/ Type-approved treatment system and its intended purpose..	42
Indications/contraindications/processing instructions.....	43
Assembly and adaption regulations	44
Labeling of the custom-made device	45
Perpedes EVA wedge.....	46
Instructions for use and informations on risks	47
Cleaning and care/storage/reuse.....	49
Disposal/guarantee.....	50

INTRODUCTION

Dear clients,

Please read the instructions for use carefully before using the product. Follow all instructions, especially the safety instructions.

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For the incorporation of foot orthoses into safety footwear in accordance with European standards DIN EN ISO 20344/20345 and 61340.

These production instructions apply for the individual incorporation of orthopedic ELTEN ESD Pro foot orthoses in ELTEN work safety shoes in accordance with compliance with DGUV Rule 112-191 (BGR 191). The individual work steps are defined and the materials and components to be used for the foot orthoses are specified here.

The instructions are required in order to ensure conformity with the EU type examination certificate for occupational safety footwear even after insertion of the foot orthoses.

Deviation from these production instructions shall render the EU type examination certificate null and void!

TYPE-APPROVED TREATMENT SYSTEM AND ITS INTENDED PURPOSE

The orthopedic ELTEN ESD Pro foot orthoses were tested and certified in combination with ELTEN work safety shoes in accordance with DIN EN ISO 20344/20345 and 61340 by recognized test institutes.

The combination of the orthopedic ELTEN ESD Pro foot orthoses and ELTEN work safety shoes is therefore type-approved and can be used in static and, depending on the shoe model, even in ESD areas.

The testing applies only for the relevant shoe models from ELTEN. The currently approved models of ELTEN work safety shoes with the orthopedic ELTEN ESD Pro foot orthoses can be found at:

www.perpedes.de/downloads/elten

INDICATIONS

The ELTEN ESD Pro foot orthoses are suitable for the conservative treatment of foot, knee, hip and lower back problems. Alongside the treatment of foot problems, they can also be used for the prevention of overloading of the entire postural and musculoskeletal system (knee, hip and back area). The treating physician is responsible for the exact diagnosis.

CONTRAINDICATIONS

- Diabetic foot syndrome

PROCESSING INSTRUCTIONS

The production instructions apply for the ELTEN ESD Pro foot orthoses with the following item numbers:

Foot orthose/Item no. GA02

ASSEMBLY AND ADAPTATION REGULATIONS

The foot orthoses may only be modified or adapted in accordance with the procedures set out below and as defined.

You can build up and adjust the inlay to meet the requirements of the wearer as long as you take account of the area beneath the toe cap which cannot be altered (the total thickness of the foot orthoses must not be changed in the area of the steel toe cap):

- Up to 10 mm shortening compensation, with or without stiffening, inner edge elevation, outer edge elevation, lowering and padding, as well as heel elevation up to a maximum of 10 mm.
- All materials required to change the foot orthoses can be freely selected if they are used below the top cover material and outside of the steel toe cap area.

LABELING OF THE CUSTOM-MADE DEVICE

The orthopedic change must be marked in the shoe below the orthopedic foot orthoses by the specialist retailer:

- a) with the name of the manufacturer
(name of the orthopedic expert)
- b) with the type number

The manufacturer of the custom-made device must apply the procedure for custom-made devices according to Annex XIII of the Regulation (EU) 2017/745 on medical devices (MDR).

PRODUCTION INSTRUCTION HAIX® CO SYSTEM

For the incorporation of foot orthoses into safety footwear in accordance with European standards DIN EN ISO 20344/20345 and 61340.

These production instructions apply for the individual incorporation of orthopedic HAIX® CO system foot orthoses in HAIX® work safety shoes in accordance with compliance with DGUV Rule 112-191 (BGR 191). The individual work steps are defined and the materials and components to be used for the foot orthoses are specified here.

The instructions are required in order to ensure conformity with the EU type examination certificate for occupational safety footwear even after insertion of the foot orthoses.

Deviation from these production instructions shall render the EU type examination certificate null and void!

TYPE-APPROVED TREATMENT SYSTEM AND ITS INTENDED PURPOSE

The orthopedic HAIX® CO system foot orthoses were tested and certified in combination with HAIX® work safety shoes in accordance with DIN EN ISO 20344/20345 and 61340 by recognized test institutes.

The combination of the orthopedic HAIX® CO system foot orthoses and HAIX® work safety shoes is therefore type-approved and can be used in static and, depending on the shoe model, even in ESD areas.

The testing applies only for the relevant shoe models from HAIX®. The currently approved models of HAIX® work safety shoes with the orthopedic HAIX® CO system foot orthoses can be found at:

www.perpedes.de/downloads/haix®

For orthopedic technicians / orthopedic shoe technicians and other qualified specialists, this results in the advantage of being able to supply in compliance with the requirements of the EU guidelines without any further application for testing procedures (subject to possible variations owing to the legal requirements in the specific country of usage).

INDICATIONS

The HAIX® CO system foot orthoses are suitable for the conservative treatment of foot, knee, hip and lower back problems. Alongside the treatment of foot problems, they can also be used for the prevention of overloading of the entire postural and musculoskeletal system (knee, hip and back area). The treating physician is responsible for the exact diagnosis.

CONTRAINDICATIONS

- Diabetic foot syndrome

PROCESSING INSTRUCTIONS

The production instructions apply for the HAIX® CO system with the following item numbers:

Foot orthose/Item no. GA01

ASSEMBLY AND ADAPTATION REGULATIONS

1. The foot orthoses may only be modified or adapted in accordance with the procedures set out below and as defined.
2. You can build up and adjust the inlay to meet the requirements of the wearer as long as you take account of the area beneath the toe cap which cannot be altered (up to 10 mm behind the edge of the toe cap - the total thickness of the foot orthoses must not be changed in the area of the steel toe cap).
3. Only an EVA wedge with a Shore hardness of at least 40° Shore A and a maximum of 50° Shore A may be used to adapt the foot orthoses. The damping wedges do not have to come from a specific manufacturer.

WARNING:

The antistatic properties and ESD capability of all EG prototype-tested HAIX® professional and safety footwear are secured via an electrically conductible textile strip on the inside of the shoe.

This grey anti-static band must extend at least 5 centimetres into the shoe. In order to continue to be able to guarantee conductivity, the inlay must not cover the strip and the contact area to the foot needs to be about 2 centimetres. In light of the above information regarding the conductive strip, the use of antistatic or ESD-capable materials (such as EVA / adhesive) is not absolutely necessary.

LABELING OF THE CUSTOM-MADE DEVICE

The orthopedic change must be marked in the shoe below the orthopedic foot orthoses by the specialist retailer:

- a) with the name of the manufacturer
(name of the orthopedic expert)
- b) with the type number

The manufacturer of the custom-made device must apply the procedure for custom-made devices according to Annex XIII of the Regulation (EU) 2017/745 on medical devices (MDR).

PRODUCTION INSTRUCTION ATLAS KLIMA-KOMFORT®

For the incorporation of foot orthoses into safety footwear in accordance with European standards DIN EN ISO 20344/20345 and 61340.

These production instructions apply for the individual incorporation of orthopedic ATLAS KLIMA-KOMFORT® foot orthoses in ATLAS work safety shoes in accordance with compliance with DGUV Rule 112-191 (BGR 191). The individual work steps are defined and the materials and components to be used for the foot orthoses are specified here.

The instructions are required in order to ensure conformity with the EU type examination certificate for occupational safety footwear even after insertion of the foot orthoses.

Deviation from these production instructions shall render the EU type examination certificate null and void!

TYPE-APPROVED TREATMENT SYSTEM AND ITS INTENDED PURPOSE

The orthopedic ATLAS KLIMA-KOMFORT® foot orthoses were tested and certified in combination with ATLAS work safety shoes in accordance with DIN EN ISO 20344/20345 and 61340 by recognized test institutes.

The combination of the orthopedic ATLAS KLIMA-KOMFORT® foot orthoses and selected ATLAS work safety shoes is therefore type-approved and can be used in static and, depending on the shoe model, even in ESD areas.

The approval only applies to atlas work safety shoes that are marked with the "insole symbol" on the Atlas homepage and in the Atlas documents. If you have any questions, please contact the manufacturer directly!

For orthopedic technicians / orthopedic shoe technicians and other qualified specialists, this results in the advantage of being able to supply in compliance with the requirements of the EU guidelines without any further application for testing procedures (subject to possible variations owing to the legal requirements in the specific country of usage).

INDICATIONS

The ATLAS KLIMA-KOMFORT® foot orthoses are suitable for the conservative treatment of foot, knee, hip and lower back problems. Alongside the treatment of foot problems, they can also be used for the prevention of overloading of the entire postural and musculoskeletal system (knee, hip and back area). The treating physician is responsible for the exact diagnosis.

CONTRAINDICATIONS

- Diabetic foot syndrome

PROCESSING INSTRUCTIONS

The production instructions apply for the ATLAS KLIMA-KOMFORT® with the following item numbers:

Foot orthose/Item no. GA04

ASSEMBLY AND ADAPTATION REGULATIONS

1. When making orthopedic foot orthoses, please use the Klima Komfort® insoles as the basis for your substructure. Only these foot orthoses were tested together with the shoe models certified according to DGUV Rule 112-191.
2. The orthopedic foot orthoses may only be placed 5 mm behind the opening of the protective toe cap and only below the insole. This is absolutely necessary to maintain the anti-static properties and for the remaining height of the toe cap. In the heel area, the total height of the orthopedic foot orthoses must not exceed 13 mm.
3. Please keep in mind that hard materials can negatively affect energy absorption. Therefore, your Material selection should not be harder than standard orthopedic cork with a Shore A hardness of 60 to 65.
4. The following adjustments are possible: shortening compensation up to 13 mm / inner and outer edge increases /lowering and/or upholstery

LABELING OF THE CUSTOM-MADE DEVICE

The orthopedic change must be marked in the shoe below the orthopedic foot orthoses by the specialist retailer:

- a) with the name of the manufacturer
(name of the orthopedic expert)
- b) with the type number

The manufacturer of the custom-made device must apply the procedure for custom-made devices according to Annex XIII of the Regulation (EU) 2017/745 on medical devices (MDR).

PRODUCTION INSTRUCTION STEITZ SECURA ERGO-SOFT ESD

For the incorporation of foot orthoses into safety footwear in accordance with European standards DIN EN ISO 20344/20345 and 61340.

These production instructions apply for the individual incorporation of orthopedic Steitz Secura ERGO-SOFT ESD foot orthoses in Steitz Secura work safety shoes in accordance with compliance with DGUV Rule 112-191 (BGR 191). The individual work steps are defined and the materials and components to be used for the foot orthoses are specified here.

The instructions are required in order to ensure conformity with the EU type examination certificate for occupational safety footwear even after insertion of the foot orthoses.

Deviation from these production instructions shall render the EU type examination certificate null and void!

TYPE-APPROVED TREATMENT SYSTEM AND ITS INTENDED PURPOSE

The orthopedic Steitz Secura ERGO-SOFT ESD foot orthoses were tested and certified in combination with Steitz Secura work safety shoes in accordance with DIN EN ISO 20344/20345 and 61340 by recognized test institutes.

The combination of the orthopedic Steitz Secura ERGO-SOFT ESD foot orthoses and selected Steitz Secura work safety shoes is therefore type-approved and can be used in static and, depending on the shoe model, even in ESD areas.

The testing applies only for the relevant shoe models from Steitz Secura. The currently approved models of Steitz Secura work safety shoes with the orthopedic Steitz Secura ERGO-SOFT ESD foot orthoses can be found at:

www.perpedes.de/downloads/steitz_secura

For orthopedic technicians / orthopedic shoe technicians and other qualified specialists, this results in the advantage of being able to supply in compliance with the requirements of the EU guidelines without any further application for testing procedures (subject to possible variations owing to the legal requirements in the specific country of usage).

INDICATIONS

The Steitz Secura ERGO-SOFT ESD foot orthoses are suitable for the conservative treatment of foot, knee, hip and lower back problems. Alongside the treatment of foot problems, they can also be used for the prevention of overloading of the entire postural and musculoskeletal system (knee, hip and back area). The treating physician is responsible for the exact diagnosis.

CONTRAINDICATIONS

- Diabetic foot syndrome

PROCESSING INSTRUCTIONS

The production instructions apply for the Steitz Secura ERGO-SOFT ESD with the following item numbers:

Foot orthose/Item no. GA03

ASSEMBLY AND ADAPTATION REGULATIONS

1. Changes in the forefoot area are possible up to 10 mm in front of the toe cap if the production instructions are observed. In the area of the toe cap (about 50 to 60 mm from the tip of the insole) no changes may be made due to the required residual height when the toe cap is loaded. The transitions must always be well sanded.
2. The required construction in the forefoot area must be made with an anti-static material and conductive adhesive that meet the requirements.
3. In addition, there is the possibility that the insole can be built up individually in the ball and joint area. The structure in the forefoot and heel area must be adhered to as described above. For better perspiration to the insole, the insole should be perforated with ERGO-SOFT ESD after completion in the forefoot area. (Hole diameter 3 mm.)
4. In the heel area, a construction height of up to 13 mm with Poro or EVA material (40 ± 5 Shore A) is possible.

LABELING OF THE CUSTOM-MADE DEVICE

The orthopedic change must be marked in the shoe below the orthopedic foot orthoses by the specialist retailer:

- a) with the name of the manufacturer
(name of the orthopedic expert)
- b) with the type number

The manufacturer of the custom-made device must apply the procedure for custom-made devices according to Annex XIII of the Regulation (EU) 2017/745 on medical devices (MDR).

APPLIES TO ELTEN ESD PRO & HAIX® CO SYSTEM

& ATLAS KLIMA-KOMFORT®

& STEITZ SECURA ERGO-SOFT ESD

PERPEDES EVA WEDGE

EVA-Worker EVA Wedge 50°Shore A, thickness available in 6.0 and 10.0 mm

ENGLISH

	6 MM EVA WEDGE (CONDUCTIVE)	10 MM EVA WEDGE (CONDUCTIVE)
Size 1 (35 bis 39)	WOK61	WOK101
Size 2 (40 bis 44)	WOK62	WOK102
Size 3 (45 bis 48)	WOK63	WOK103
Antistat adhesive Ultraflex, 1 kg	79AK1	79AK1
Thinner for Antistat, 1 liter	79IS1	79IS1

INSTRUCTIONS FOR USE AND INFORMATION ON RISKS

- The prefabricated foot orthoses must be properly customized for the foot and adjusted in the shoe by a qualified specialist.
- Foot orthoses should be worn with (crease-free) socks and should always be used as a pair! Exceptions must be stated by the treating physician.
- The shoes must ensure a sufficient interior volume (shoe length / width, toe space) for the foot as well as for the chosen prefabricated foot orthoses. The inlays included in the scope of delivery for the occupational safety footwear must be removed before the prefabricated foot orthose is inserted.
- The chosen prefabricated foot orthoses need to be cut to match the width and length. The forefoot, mid foot and rear foot must have sufficient flat support from the foot orthoses and should not be on the edge as this can lead to the formation of blisters or calluses on the soles of the feet.
- Clear identification of the left and right foot orthoses must be ensured for the patient / user.
- The foot orthoses should not be exposed to temperatures above 60 °C (e.g. heating, driers, direct sunlight).
- The foot orthoses are not suitable for direct contact with acids and bases. The manufacturer provides no warranty in this case.
- Materials which are used for the foot orthoses are subject to aging and their properties may change. If the coverings, cushioning material or orthotic shells have become damaged, holey or wrinkled as a result of mechanical stress or are broken then the foot orthoses should no longer be worn!

Foot orthoses are affected by body weight, pressure, torsion and rubbing and therefore have a load-dependent lifespan. Temperatures, sock materials, orthotic materials, footwear and environmental factors can also limit the functional life.

- Checking of the function and condition of the foot orthoses by an orthopedic technician / orthopedic shoe technician or other qualified specialist after a maximum of 6 months is recommended.

CLEANING AND CARE

- Foot orthoses must be removed from the shoes daily in order for moisture to evaporate.
- The surface materials can be cleaned as required with mild washing lotion and a damp sponge or cloth. Soap residues should be removed. Do not use any external heat sources to dry, simply leave to air-dry. Note that the foot orthoses are not machine-washable!
- Mild disinfectant (suitable for human skin) can be sprayed on / used in small amounts. Discoloration of the covering material as a result cannot be excluded and is not grounds for a warranty claim.

STORAGE

- Store in the sealed original packaging between 10 °C and 25 °C without mechanical influences and avoiding direct UV radiation.
- Regular checking of the storage and conditions is recommended.

REUSE

- The product is used for an individually fitted treatment for an individual foot orthoses wearer and should not be reused or subsequently used by other third persons.
- The manufacturer accepts no liability in the event of improper use.

DISPOSAL

You may dispose of the product in accordance with the applicable national regulations.

GUARANTEE

The current statutory provisions apply. Disregard of the instructions for use, improperly performed cleaning as well as independent technical changes and additions to the product will void the warranty and the general product liability.

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Perpedes GmbH